



# LM-1200 / 1000(P)

AUTO LENSMETER



Eye & Health Care

# LM-1200 / 1000(P) Auto Lensmeters


*Advanced & Dynamic Technologies Provide for Faster, Easier and More Accurate Measurement*

## ***Faster & Easier for Progressive Power Lenses***

### **HARTMANN SENSOR with 108 Multiple Measuring Points**


Advanced simultaneous measurement of 108 multiple points within the nosepiece provides easier and faster measurement with greater accuracy and reliability.

### **HARTMANN SENSOR**




### **Conventional**

Unable to determine the direction of the reading point without moving the lens around.




Unable to detect the reading point immediately.




### **LM-1200 / 1000(P)**

Instantly decides the direction of the reading point.



Detects the reading point immediately.



# LM-1200 / 1000(P) Auto Lensmeters

## *Enhanced Functionality, Operability & Versatility*

### ■ **World's First Measurement of High Power Progressive Lenses**


The prism measurement range is expanded to 20Δ, offering greater versatility.

### ■ **Automatic Lens Type Detection**


Just place the lens on the nosepiece, and the system detects the lens type - single vision lens or progressive power lens - and switches its measuring mode automatically.

### ■ **UV Transmittance Measurement**

UV transmittance of the lens is graphically displayed for 4 levels. The measured data can be saved, and the UV transmittance of multiple lenses can be compared.



Comparison



### ■ **Improved Marking Dots**

Newly developed marking ink makes it easy to mark and provides clear dots even on lenses with water repellent coating / finish.



# LM-1200 / 1000(P) Auto Lensmeters

## Ergonomic & User-friendly Design

### Ergonomic Design

The new user-friendly ergonomic design provides smooth measurement.



Easy measurement even for glasses without hinges.

### High-Speed Line Printer (LM-1200 / 1000P)

The LM-1200 and LM-1000P provide fast and easy-to-read printouts. A model without the printer module (LM-1000) is also available.

#### Sample Printout (LM-1200)

〈PROGRESSIVE〉			
RIGHT		LEFT	
+1.50	SPH	+1.25	
-0.25	CYL	+0.00	
137°	AXS	0°	
1.50	PSM	3.25	
268°	BAS	266°	
+2.75	ADD	+2.75	
17	LEN	17	Progressive length*
LEVEL1	UV	LEVEL1	UV transmittance
----- WIDTH / LEN -----			
20 / 17		18 / 17	Measurement position from channel width / distance portion*
----- PD -----			
32.0	63.5	31.5	Right PD / Binocular PD / Left PD*
1.5	INS	-2.0	Near portion inside amount*
NIDEK LM-1200			

\*LM-1200 only

### User-Friendly Tilttable LCD

The tilttable (30°) full-graphic LCD monitor provides easier operation for both standing and sitting operators.

### Small Footprint

The LM-1200 / 1000(P) is compact and space saving.

### USB Interface

Provides for data upload to PCs and servers.



# LM-1200 Outstanding Functions


## Easy PD Measurement

The LM-1200 can measure the pupillary distance (the distance between the optical centers) easily - for both near and distance portions.



Offers automatic Right / Left detection with the special PD slider, which can also help the operator easily hold glasses while measuring.

## Addition Power Change Graph



The LM-1200 can also graph the power transition of progressive lenses to indicate if the power remains the same, increases, or decreases around the near portion. This function is particularly useful when it is necessary to duplicate the customer's old lenses.

## Color LCD Monitor

The LM-1200 provides a clear color LCD monitor.



## LM-1200 / 1000(P) Specifications

Measurement Range	LM-1200	LM-1000 / 1000P
Sphere (Spectacle lenses) (Contact lenses)	-25.0 D to +25.0 D -25.0 D to +25.0 D (BC=6.00 to 9.00) (0.01 / 0.06 / 0.12 / 0.25 D increments)	
Cylinder	0 D to $\pm 9.99$ D (-, MIX, +) (0.01 / 0.06 / 0.12 / 0.25 D increments)	←
Axis	0° to 180° (1° increments)	
ADD	0 D to $\pm 9.99$ D (Add and Ad2) (0.01 / 0.06 / 0.12 / 0.25 D increments)	
Prism	0 $\Delta$ to 17 $\Delta$ (horizontal), 0 $\Delta$ to 20 $\Delta$ (vertical), (0.01 / 0.06 / 0.12 / 0.25 $\Delta$ increments)	
Prism mode	$\Delta$ $\theta$ , BI/O BU/D	←
PD Measurement	40 - 100 mm (0.5 mm increments)	None
Measurable Lens diameter	$\phi$ 5 to 120 mm	←
Measuring time	0.13 sec. (Minimum)	←
Measurable transmittance	10% and over (20% and over for $\pm 15.0$ D to $\pm 25.0$ D)	←
UV transmittance	4 levels (None, Low, Moderate, High) with 365 nm (UV-A)	←
Wavelength / Measuring point	660 nm (Red) / 108 within nosepiece	←
Display	30° tilt Color LCD with back light	30° tilt LCD with back light (Black / White)
Printer	High speed built-in line printer (Paper width: 58 mm)	← (LM-1000P only)
Interface	RS-232C, USB	←
Marking system	Ink cartridge type (or ink pad type)	←
Power source	AC100 V to 120 V / 200 V to 240 V 50 / 60 Hz	←
Power consumption	40 VA	←
Dimensions / Weight	213 (W) x 428 (H) x 227 (D) mm / 7.5 kg 8.4 (W) x 16.85 (H) x 8.94 (D)" / 16.6 lbs.	213 (W) x 428 (H) x 227 (D) mm / 7.2 kg 8.4 (W) x 16.85 (H) x 8.94 (D)" / 15.9 lbs.
Standard accessories	Printer paper (x3), Dust cover (x1), Fuse (x2), Power cord (x1), Contact lens nosepiece (x1)	←
Optional accessories	RS-232C interface cable, USB cable (with driver), foot switch, marking cartridge (Blue & Red), EyeCare card system	←


**NIDEK CO., LTD.**
**HEAD OFFICE**

34-14 Maehama, Hiroishi  
Gamagori, Aichi 443-0038, Japan  
Telephone : 81-533-67-6611  
Facsimile : 81-533-67-6610  
URL : <http://www.nidek.co.jp>

**TOKYO OFFICE**

(International Div.)  
6F Takahashi Bldg.,  
3-2 Kanda-Jinboucho  
Chiyoda, Tokyo 101-0051, Japan  
Telephone : 81-3-3288-0571  
Facsimile : 81-3-3288-0570  
URL : <http://www.nidek.com>

**NIDEK INC.**

47651 Westinghouse Drive  
Fremont, CA 94539, U.S.A.  
Telephone : 1-510-226-5700  
: 1-800-223-9044 (US only)  
Facsimile : 1-510-226-5750  
URL : <http://www.usa.nidek.com>

**NIDEK TECHNOLOGIES AMERICA INC.**

5500 West Friendly Ave.  
Suite 101  
Greensboro, NC 27410, U.S.A.  
Telephone : 1-336-851-0225  
: 1-888-382-5064 (US only)  
Facsimile : 1-336-851-0917  
URL : <http://www.nidektech.com>

**NIDEK SOCIÉTÉ ANONYME**

Europarc  
13, rue Auguste Perret  
94042 Créteil, France  
Telephone : 33-1-49 80 97 97  
Facsimile : 33-1-49 80 32 08  
URL : <http://www.nidek.fr>

**NIDEK TECHNOLOGIES SRL**

Via dell'Artigianato, 6 / A  
35020 Albignasego (Padova), Italy  
Telephone : 39 049 8629200 / 8626399  
Facsimile : 39 049 8626824  
URL : <http://www.nidektechnologies.it>

**VISIONARY PERFORMANCE**

\*Specifications and design are subject to change without notice for improvement.



Printed on environment-friendly recycled paper.