



REFRACTOR

RT-600



Eye & Health Care

Operation, Accuracy --- Features that have gained the trust of our customers



Smooth, accurate measurement of spherical and cylindrical features

Measurement of spherical features can be made in the range between +16.75 and -19.00 D; for cylindrical features, between 0 and -6.00 D, each adjustable in increments of 0.25 D. Even greater accuracy can be obtained by making measurements in 0.12 D increments by using the auxiliary lenses, which permit an increased range of between +26.75 and -29.00 D (Option) for spherical measurements, and up to -8.00 for cylindrical measurements. In addition, the spherical power quick forwarding dial ① can be used to make adjustments in 3.00 D step intervals, resulting in efficient, as well as accurate measurements.

Unique cross cylinder configuration

Both the cylindrical axis and cross cylinder lens can be engaged, thanks to the unique NIDEK configuration, which permits easy control with a simple turn of the cross cylinder knob ②. Three types of ± 0.25 D, ± 0.37 D, and 0.50 D cross cylinder lens are available. Both ± 0.37 D and ± 0.50 D are options.

Auxiliary lenses

A full complement of auxiliary lenses are built-in ③, which are enable to do any type of measurement.

Convergence system

Shifting the convergence lever ④ to the inside lets the patient be examined under optimal conditions.

Specifications

Measurement Range

Sphere	+16.75 - -19.00 D When using aux. Lens (+/-10 D: Option): +26.75 - -29.00 D
Cylinder	0.00 D - -6.00 D When using aux.lens (-2.00 D): 0.00 D - -8.00 D
Axis	0° - 180°
Prism	0 - 20 Δ

Measurement Increments

Sphere	0.25 D When using aux. Lens (+0.12 D): 0.12D
Cylinder	0.25 D When using aux. Lens (-12.0 D): 0.12 D
Axis	5°
Prism	1 Δ

Cross Cylinder Configuration

Operates cylindrical axis
Cross cylinder lens: +0.25,
 ± 0.37 D (Option), ± 0.50 D (Option)

PD Adjustment

Convergence Adjustment

Forehead Rest Adjustment

50 - 75 mm
 ∞ -380 mm (when PD is 64 mm)
16 mm

Vertex Distance

Dimensions & Weight

13.75 mm (Normal)
336 (W) x 94 (D) x 293 (H) mm
Approx. 5 kg

Auxiliary Lenses

	Symbol	Name (details)	
For right eye	O	Open aperture	
	R	Retinoscope lens +1.50 D (standard)	
	P	Polaroid (polarization sheet) 135°	
	RMV	Vertical Maddox rod (red)	
	RMH	Horizontal Maddox rod (red)	
	RL	Red lens	
	⊙	Test mark for PD adjustment	
	+12	Sphere lens +0.12 D	
	PH	Pin hole	
	6 Δ U	6 Δ base up	
	$\pm .50$	$\pm .50$ D fixed cross cylinder	
	OC	Occluder	
	For left eye	O	Open aperture
		R	Retinoscope lens +1.50 D (standard)
P		Polaroid (polarization sheet) 45°	
WMV		Vertical Maddox rod (white)	
WMH		Horizontal Maddox rod (white)	
GL		Green lens	
⊙		Test mark for PD adjustment	
+12		Sphere lens +0.12 D	
PH		Pin hole	
10 Δ I		10 Δ base in	
$\pm .50$ D		$\pm .50$ D fixed cross cylinder	
OC		Occluder	

Specifications and design are subject to change without notice for improvement.



HEAD OFFICE

:34-14, Maehama, Hiroishi Gamagori, Aichi 443-0038, Japan
Telephone :81-533-67-6611
Facsimile :81-533-67-6610
URL :http://www.nidek.co.jp

NIDEK SOCIÉTÉ ANONYME

:Europarc 13, rue Auguste Perret 94042 Créteil, France
Telephone :33-1-49 80 97 97
Facsimile :33-1-49 80 32 08
URL :http://www.nidek.fr

TOKYO OFFICE (International Div.)

:6F Takahashi Bldg., 3-2 Kanda-Jinboucho Chiyoda,
Tokyo 101-0051, Japan
Telephone :81-3-3288-0571
Facsimile :81-3-3288-0570
URL :http://www.nidek.com

NIDEK TECHNOLOGIES SRL.

:Via dell'Artigianato, 6/A 35020 Albinasego (Padova), Italy
Telephone :39 049 8629200 / 8626399
Facsimile :39 049 8626824
URL :http://www.nidektechnologies.it

NIDEK INC.

:47651 Washington Drive Fremont, CA 94539, U.S.A.
Telephone :1-510-226-5700
:1-800-223-9044 (US only)
Facsimile :1-510-226-5750
URL :http://www.usa.nidek.com

NIDEK TECHNOLOGIES AMERICA INC.

:5500 West Friendly Ave, Suite 101 Greensboro, NC 27410, U.S.A.
Telephone :1-336-851-0225
:1-888-382-5064 (US only)
Facsimile :1-336-851-0917
URL :http://www.nidektech.com